



Installation and Operation Manual

Document Number 000-vivo-110-02-201208

Offices: **Ocean Optics, Inc. World Headquarters**
830 Douglas Ave., Dunedin, FL, USA 34698
Phone 727.733.2447
Fax 727.733.3962
8 a.m.– 8 p.m. (Mon-Thu), 8 a.m.– 6 p.m. (Fri) EST

E-mail: **Info@OceanOptics.com** (General sales inquiries)
Orders@OceanOptics.com (Questions about orders)
TechSupport@OceanOptics.com (Technical support)



WARNING

Protective Eye Wear Must Be Worn When
Using This Instrument -
Intense Ultraviolet Radiation Present
See Important Safety Notices inside.

—A—
**HALMA
GROUP
COMPANY**

**Additional
Offices:**

Ocean Optics Germany GmbH

Maybachstraße 11, D-73760

Ostfildern, Germany

Phone +49 (0)711 341696-0

Fax +49 (0)711 341696-85

E-Mail Info@OceanOptics.eu

Ocean Optics Asia

666 Gubei Road, Kirin Tower, Suite 601B, Changning District,
Shanghai, PRC. 200336

Phone 86.21.5206.8686

Fax 86.21.5206.8686

E-Mail Sun.Ling@OceanOptics.com

Ocean Optics EMEA

Sales and Support Center

Geograaf 24, 6921 EW DUIVEN, The Netherlands

Phone 31-26-3190500

Fax 31-26-3190505

E-Mail Info@OceanOptics.eu

Copyright © 2012 Ocean Optics, Inc.

All rights reserved. No part of this publication may be reproduced, stored in a retrieval system, or transmitted, by any means, electronic, mechanical, photocopying, recording, or otherwise, without written permission from Ocean Optics, Inc.

This manual is sold as part of an order and subject to the condition that it shall not, by way of trade or otherwise, be lent, re-sold, hired out or otherwise circulated without the prior consent of Ocean Optics, Inc. in any form of binding or cover other than that in which it is published.

Trademarks

Microsoft, Windows, Windows 95, Windows 98, Windows Me, Windows NT, Windows 2000, Windows XP and Excel are either registered trademarks or trademarks of Microsoft Corporation.

Limit of Liability

Every effort has been made to make this manual as complete and as accurate as possible, but no warranty or fitness is implied. The information provided is on an "as is" basis. Ocean Optics, Inc. shall have neither liability nor responsibility to any person or entity with respect to any loss or damages arising from the information contained in this manual.

Important Safety Notices

1. Read this manual before you attempt to use this instrument.
2. Do not remove or modify any installed safety device on this equipment. Doing so will void your warranty and create an unsafe operating environment.
3. Only allow qualified personnel to service this unit.
4. Before using the optical tool for the first time check for transport damage. Do not use the unit if it is damaged in any way. Contact your dealer for repair or replacement information.

Table of Contents

| | |
|---|----------|
| About This Manual | iii |
| Document Purpose and Intended Audience..... | iii |
| Document Summary..... | iii |
| Product-Related Documentation | iii |
| Upgrades..... | iii |
| Chapter 1: Setup and Operation | 1 |
| Overview..... | 1 |
| Unpacking Vivo | 2 |
| Package Contents | 2 |
| Vivo Lamp Controls..... | 2 |
| Connecting an Optical Fiber | 3 |
| Changing the Bulbs..... | 4 |
| Chapter 2: Vivo Specifications | 7 |
| Vivo Parts | 7 |
| Index | 9 |

About This Manual

Document Purpose and Intended Audience

This document provides you with an installation section to get your light source up and running.

Document Summary

| Chapter | Description |
|--|---|
| Chapter 1: Setup and Operation | Provides a list of package contents and unpacking instructions. Also contains fiber connection and bulb replacement instructions. |
| Chapter 2: Vivo Specifications | Contains product specifications and a parts list. |

Product-Related Documentation

You can access documentation for Ocean Optics products by visiting our website at <http://www.oceanoptics.com>. Select *Technical* → *Operating Instructions*, then choose the appropriate document from the available drop-down lists. Or, use the **Search by Model Number** field at the bottom of the web page.

You can also access operating instructions for Ocean Optics products on the *Software and Technical Resources* CD included with the system.

Engineering-level documentation is located on our website at *Technical* → *Engineering Docs*.

Upgrades

Occasionally, you may find that you need Ocean Optics to make a change or an upgrade to your system. To facilitate these changes, you must first contact Customer Support and obtain a Return Merchandise Authorization (RMA) number. Please contact an Ocean Optics Application Scientist for specific instructions when returning a product.

Setup and Operation

Overview

The Vivo is an NIR light source for analysis of pharmaceutical products, food and more. Vivo features four tungsten-halogen bulbs placed at a 90 degree angle from your optical fiber that can be turned on and off individually. By providing output from four tungsten-halogen bulbs, Vivo becomes a high-powered VIS-NIR light source, enabling you to integrate your spectrometer for a shorter time period than conventional light sources.

You can attach Vivo directly to an Ocean Optics Stage-RTL-T, or to your own stand for great stability and control. Also, Vivo features active cooling that reduces the risk of overheating your sample, ensuring accurate measurements.



Vivo

Unpacking Vivo

The following sections provide instructions on unpacking and setting up your Vivo.

Before using the Vivo for the first time, check for transport damage. Be sure to adhere to all warnings on the unit and in this manual.

► Procedure

1. Unpack your Vivo device carefully.
2. Inspect the outside of the instrument and make sure that there is no damage. Do not use the instrument if damage is present.
3. Switch on the unit and check if all bulbs are working by pressing both pushbuttons.

WARNING

Do not stare directly into the light beam of the lamps!

Package Contents

Your Vivo package should contain the following:

- ❑ Vivo device
- ❑ Power supply (12 Vdc @ 2.5 A)
- ❑ Three hexagonal socket screw keys (SW 1.3 mm, SW 2.0 mm and SW 2.5 mm)
- ❑ SMA barrel
- ❑ Magnetic ring

Caution

For optimal use, take care that no objects or liquids enter the device and that there is no grease residue on the glass window.

Vivo Lamp Controls

The device is running with 12 Vdc @ 2.5 A maximum current. Be sure to use the included power supply and do not change its voltage settings. Other voltages may cause damage.

At the right side of the device there is a green power switch next to the DC power socket. On the front there are two push buttons and three rotating switches. The push buttons can be used to switch on one pair of bulbs. The rotating switches regulate the intensity of the bulbs. Refer to the picture below.




Connecting an Optical Fiber

To connect an optical fiber with an SMA connection, use the SMA-barrel. The fiber with the SMA-barrel must be mounted from the bottom of the device.

► Procedure

1. Push the fiber till the SMA-barrel completes with the end of the SMA lead-through.



2. To fasten the SMA-barrel, tighten the screw shown below with a hexagonal socket screw key SW 2.5.



Changing the Bulbs

Caution

Bulbs must be changed as a pair (two face-to-face bulbs are considered to be a pair).

Do not over-tighten the screws. The plastic screws could be damaged and the stud screws could damage the bulbs if they are tightened too much.

► Procedure

1. Before changing the bulbs switch off the power and disconnect the power plug.
2. Unfasten the lamps by dismounting the top plate. To do this, unscrew the 4 plastic screws on the top plate (shown in figure below).



3. Disassemble the punched plate on the bottom side by removing the 4 screws with the hexagonal socket screw key SW 2.0.
4. After removing the punched plate, unplug the lamps you want to change. The lamp lines are yellow and black, and the plugs are placed in the corners of the circuit board (see figure below).



5. To remove a bulb, unfasten a stud screw from the top side (see figure in Step 2) using a hex screw key SW 1.3.
6. Pull the bulb on its line to the button or push it to the top to remove the bulb. Replace the new bulb in the same way. The position of the new bulb should be the same as shown in the figures below. Do not place the bulbs nearer to the glass window (this could cause stray light or cover the fiber) or farther (this would defocus the light spot on the glass window and the spot becomes less uniform).



Chapter 2

Vivo Specifications

| Specification | Value |
|------------------------|-----------------------------------|
| Size (LxHxW) in mm | 120 x 160 x 75 |
| Weight | 1.6 kg |
| Wavelength Range | 360 – 1700 nm |
| Stability | 0.15 % @ 900 nm; 0.05 % @ 1400 nm |
| Drift | < 0.2 % per hour |
| Time to stabilize | 2 – 5 minutes |
| Lamp power | 30 W |
| Output to bulb | up to 5 VDC / 0.97 A |
| Bulb lifetime | 2000 hours |
| Bulb color temperature | 2800K |
| Temperature range | 5 -- 35°C |
| Humidity range | 5 – 95% at 40°C |

Vivo Parts

| Part | Catalog Number |
|----------------------------------|----------------|
| Vivo IR Illuminator Light Source | VIVO |
| Vivo replacement bulb | VIVO-B |

Index

| | | | |
|------------------|----------|------------------------------------|----------|
| | B | | O |
| bulb | | operation, 1 | |
| changing, 4 | | | P |
| | D | package contents, 2 | |
| document | | parts, 7 | |
| audience, iii | | product-related documentation, iii | |
| purpose, iii | | | S |
| summary, iii | F | safety notices, A | |
| | | setup, 1 | |
| fiber | | specifications, 7 | |
| connecting, 3 | L | | U |
| | | unpacking, 2 | |
| lamp controls, 2 | | upgrades, iii | |

